




## PF-BT Cartridge Filter Data Sheet



### Features

- Low Pressure Drop • High Flow • Long Service
- Steam Sterilization while left installing to the tank
- Compatibility with existing Vent Filter
- Excellent Bacteria Removal Performance

### Application

- For the process of Food and Beverage, Pharmaceutical company, Cosmetic Industry, Chemical Company

### Specification

#### ● Material

Membrane: PTFE End-cap•Core and Support Material: PP O-Ring: Silicone

#### ● Dimension O.D.:70mm I.D.:32mm Length: 288mm

#### ● Membrane Surface Area: 7000cm<sup>2</sup>

#### ● Pressure-Resistant : 0.35MPa (20°C)

#### ● Heat-Resistant : SIP 121°C Half-Hour 100times

#### ● Bacteria Removal Performance :

Brevundimonas diminuta (ATCC 19146)


Aerosol-Challenge-Test LRV  $\geq 9$

#### ● 100% Tested and Proved to be Integrity Test

#### ● 1pc / 1 box

#### ● Ordering Guide

BT-01T K72 32-70-288



\* Description of each characteristic data is based on our original test methods. The actual data may differ slightly from the data shown here by use conditions.

\* Please be reminded that we may change this content without advance notice.

\* Copyright JNC Filter Co., Ltd. All rights reserved. No reproduction or republication without written permission.



JNC Filter Co., Ltd.

Headquarters and Osaka Sales Office  
Tokyo Sales Office  
Nagoya Sales Office  
web : <http://www.jnc-corp.co.jp/filter/>

Tel : +81-6-6441-3257 Fax : +81-6-6441-8836  
Tel : +81-3-3243-6222 Fax : +81-3-3243-6223  
Tel : +81-52-583-6007  
Email : [jft.exp@jnc-corp.co.jp](mailto:jft.exp@jnc-corp.co.jp)

3-3-23, Nakanoshima, Kita-ku, Osaka 530-6108  
2-2-1, Otemachi, Chiyoda-ku, Tokyo 100-8105  
4-4-10, Meieki, Nakamura-ku, Nagoya 450-8561

